


bright Stent Graft

Exceptional conformability and flexibility


Multiple X-markers
for precise placement

Contoured edge
for precise graft integration

GENOSS



bright Stent Graft



Characteristics

- Exceptional Flexibility
- Reliable Sealing Conformability
- Accurate & Predictable Deployment


Conformability / Flexibility



GENOSS




GENOSS




Company A

Radial force



GENOSS	Company A
34.1 N	25.9 N

3-point bending



GENOSS	Company A
0.089 N	0.553 N

Order Information

Stent Diameter (mm)	Stent Length (mm)	Product Codes	Recommended Introducer(Fr)	Stent Diameter (mm)	Stent Length (mm)	Product Codes	Recommended Introducer(Fr)
6.0	40	GSG-06-040	8	8.0	40	GSG-08-040	8
	60	GSG-06-060			60	GSG-08-060	
	80	GSG-06-080		9.0	40	GSG-09-040	9
7.0	40	GSG-07-040			60	GSG-09-060	
	60	GSG-07-060			10.0	40	
	80	GSG-07-080					

L_Stent Graft(Rev.0,2603)